



Measurement of Magnetic Particle Concentrations in Wildfire Ash via Compact NMR

Jacob Martin


Department of Mech. Eng.

Department of Physics and Astronomy


University of South Carolina

Why monitor magnetic contents of wildfire ash?

- Effects on topsoil
 - Ash deposits enhance magnetic content in soil
 - Magnetic properties are closely related to climate & rainfall
- Deposition through runoff water
 - Nearby bodies of water accumulate magnetic content
 - Nanoscale magnetite is linked to brain disease
- Understand fire severity and the reaches of magnetic deposition



USGS, "How wildfires threaten U.S. water supplies," Water Data Labs, 06-Nov-2020. [Online]. Available: [https://labs.waterdata.usgs.gov/visualizations/fire-hydro/index.html#/.](https://labs.waterdata.usgs.gov/visualizations/fire-hydro/index.html#/) [Accessed: 28-Oct-2022].




Goals

- Quickly assess total magnetic particle (MP) concentrations in ash & water
 - Magnetite + maghemite
 - Public health school currently uses synchrotron (complex, time-consuming, not portable)
- Develop a system for eventual in situ use
- Best candidate: nuclear magnetic resonance (NMR)

Nuclear magnetic resonance (NMR) techniques

- High-field NMR spectroscopy


- High quality molecular structure data
- Information about environment around nucleus
- High cost, not portable, difficult to achieve magnetic field homogeneity



B. Blümich, "Introduction to compact NMR: A review of methods," TrAC Trends in Analytical Chemistry, vol. 83, pp. 2–11, Oct. 2016.

- Low-field NMR relaxometry



- Very low-quality molecular structure data
- Information about **environment around nucleus (relaxation rate)**
- Low cost, portable, **no demands on field homogeneity**



B. Blümich, "Introduction to compact NMR: A review of methods," TrAC Trends in Analytical Chemistry, vol. 83, pp. 2–11, Oct. 2016.

NMR Relaxometry with MPs



- Transverse relaxation (R2) rates are related to MP concentration (magnetic environment)
- Linear relationship is well established
- **Objective: Dissolve ashes in water to measure MP content via compact NMR**



Gunn, J., Paranjli, R. K., and Zhang, M., 2009. "A simple and highly sensitive method for magnetic nanoparticle quantitation using ¹H-NMR spectroscopy". *Biophys J*, 97(9), Nov., pp. 2640–2647.

Permanent Magnet Array


- Components:
 - N42 NdFeB permanent magnets (8)
 - 1018 steel return yokes
 - 1018 steel magnet caps (2)
- 2D finite element simulation
 - Field strength ≈ 0.65 T



Permanent Magnet Array

- **0.565 T** strength at 23°C
 - -800 ppm/K gradient
- Larmor (operating) frequency:
 - $f_{Larmor} = \gamma B = \left(42.58 \frac{MHz}{T}\right) (0.565 T) \approx \mathbf{24 MHz}$
- 150 ppm homogeneity
- 4.4 lbs


N42 magnet




N42 magnet



Sample
location




Fully assembled



3D printed guide

Signal Processing & Amplification Electronics



- All ports & PCB traces matched to 50 Ω
- One 24 V power supply required
- 24 MHz pulse amplified to 35 dBm
- Duplexer isolates probe and LNA
- NMR signal amplified 80 dB before acquisition





Signal Generation, Control, & Data Acquisition


- NI PXI chassis
 - Arbitrary waveform generator
 - Pulse train generator
 - 16-bit digitizer




- CPMG pulse sequence
 - 90° pulse time = 6 μ s
 - τ = 0.625 ms

Signal Generation, Control, & Data Acquisition

- Thermocouple used for frequency calibration (magnet gradient)
- LabVIEW GUI
- 8 scans (averages) takes 2 min.
- Transverse decay model:
 - $M_{xy} = M_0 \exp(-R_2 t)$




Full Desktop System




- 10 total ash samples
- 20 mg in 20 mL of water
- Distilled water used as reference
- R_2 extracted via least squares regression


- $M_{xy} = M_0 \exp(-R_2 t)$





- Ash samples analyzed using synchrotron
- Total iron content collected
- Only magnetic iron (**maghemite & magnetite**) affect R2 rates
- Linear relationship verified





- Magnetic content in wildfire ash can be measured via compact NMR
- Robust and compact permanent magnet array sufficient for 5mm NMR tubes
- Simple and rapid tests to quantify MP content compared to current method
- Currently working to implement for in situ monitoring

Compact-NMR

intensity

t (s)

University of South Carolina

Thank you!

Questions?