

Assessing Magnetic Particle Content in Algae Using Compact Time Domain Nuclear Magnetic Resonance

Parker Huggins¹, Win Janvrin¹, Jake Martin^{1,3}, Ashley Womer⁴, Austin Downey^{1,5}, John Ferry⁴,
Mohammed Baalousha⁶, Jin Yan⁷

¹Department of Mechanical Engineering, University of South Carolina


²Department of Physics and Astronomy, University of South Carolina

³Department of Chemistry and Biochemistry, University of South Carolina

⁴Department of Civil and Environmental Engineering, University of South Carolina

⁵Department of Environmental Health Sciences, University of South Carolina

⁶SRI, Palo Alto, CA



Satellite image of phytoplankton swirling around the Swedish island of Gotland in the Baltic Sea, in 2005

NASA Goddard Space Flight Center
Credit: USGS/NASA/Landsat 7 - Flickr:
Van Gogh from Space



UNIVERSITY OF
SOUTH CAROLINA



Harmful algal blooms in Harford
County, Md
Chesapeake Bay Program

Introduction

Materials &
Methods

Results &
Discussion

Conclusion

Algae blooms

- Algae plays an important role in ecosystem wellness
 - Form the base of aquatic food webs
- Algae overgrowth, however, can be detrimental to an ecosystem
 - Introduce toxins that affect the availability of safe drinking water
 - Block sunlight needed for aquatic organisms
 - Contribute to the depletion of a habitat's oxygen levels





Photo Credit: Dr. Jennifer L. Graham | U.S. Geological Survey




Why assess magnetic particle (MP) content in algae?

- Iron uptake mechanisms of algae have been the subject of much research
 - Especially pertinent in commercial applications requiring the growth of microalgae
- A quick & reliable method for monitoring MP concentrations in algae has an abundance of applications
 - Water quality monitoring
 - Conservation initiatives
 - Iron harvesting



Iron-oxidizing bacteria in surface water

NH Estuaries Project - taken by
New Hampshire Estuaries Project
(www.nhep.unh.edu)



Goals (now and future)

1. Quickly assess the MP content of algae & surrounding water
2. Monitor the magnetic iron uptake of algae blooms/mats
3. Develop a system for in situ MP estimations

Approach: Time domain nuclear magnetic resonance (TD-NMR) to monitor iron uptake in *Lyngbya Wollei* (recently renamed *Microseira wollei*)




Filamentous cyanobacterium of a genus *Lyngbya*, as collected in Baja California


NASA - <http://microbes.arc.nasa.gov/images/content/gallery/lightms/publication/lyngbya.jpg>

Nuclear magnetic resonance (NMR) techniques

- High-field NMR spectroscopy
 - Frequency domain analysis
 - High resolution
 - Expensive & bulky
- Low-field NMR relaxometry
 - Time domain analysis
 - Low resolution
 - Inexpensive and portable



B. Blümich, "Introduction to compact NMR: A review of methods," TrAC Trends in Analytical Chemistry, vol. 83, pp. 2–11, Oct. 2016.



B. Blümich, "Introduction to compact NMR: A review of methods," TrAC Trends in Analytical Chemistry, vol. 83, pp. 2–11, Oct. 2016.

Introduction


Materials &
Methods


Results &
Discussion

Conclusion

ARTS-Lab desktop NMR system



- Control handled by LabVIEW program and NI-PXI chassis
- All electronics (barring two amplifiers) housed on a single PCB
- GUI developed for easy data acquisition and export






Permanent magnet array

- N42 cylindrical dipole magnets enclosed by a steel yolk
- 1018 carbon steel caps affixed to magnet surfaces
- Peak flux density of 0.645 T \rightarrow Larmor frequency of 27.5 MHz
- Temperature shift gradient of -800 ppm/K







RF electronics

- A single 24 V DC power supply required
- Impedance of all cables and PCB traces matched to $50\ \Omega$
- Waveform generator \rightarrow sine wave at Larmor frequency
- Pulse generator \rightarrow CPMG pulse train
- Duplexer (crossed diodes) isolates probe and LNA


General flow







Signal generation and control

- NI PXI chassis
 - Arbitrary waveform generator
 - Pulse train generator
 - 16-bit digitizer




- CPMG pulse train
 - 3955 total pulses
 - 90° pulse duration is $6 \mu\text{s}$
 - $\tau = 0.625 \text{ ms}$







Data acquisition


- LabVIEW GUI serves as front end
- Each test comprises 5 scans (averages)
- Time for T_2 curve acquisition < 1 min
- Thermocouple used for frequency calibration



TD-NMR signals and MP content

- T_2 relaxation modeled as $M_{xy}(t) = M_0 \exp(-t/T_2)$
- Relaxation rate is the reciprocal of relaxation time (i.e., $R_2 = 1/T_2$)
- Linear relationship between R_2 and MP concentration well established





Sample collection

- Four algae samples were sourced from Lake Wateree (South Carolina)
- Samples collected on different dates
- Three samples (A-01, A-02, A-03) frozen with liquid nitrogen and submerged in MQ water
- Fourth sample (A-04) mixed directly with MQ water and sonicated




(a)




(b)

Relaxation rates

- Various locations along NMR tubes housing algae samples were tested
- Large relaxation rates observed (comparable to those of magnetite samples)
- Decay rates maximized in algae-dense regions






MP separation


- Fourth algae sample tested before & after undergoing MP separation
- MP separator comprises N52 permanent magnet and 3D printed housing
- Large decrease in decay rate observed following MP separation



(a)





(b)



Preliminary water data

- Four water samples (algae free) gathered at increasing distances from an algae mat
- Samples tested 5 times with each test comprising 5 individual scans
- T_2 time increases with distance from the algae mat






Algae and MPs

- Compact NMR shows promise for MP-based estimations concerning algae
- Clear relationship between R_2 and MPs observed
- Proposed scheme can rapidly assess the relative MP concentration of samples

Future work

- Flow-through system for accelerated probing
- Electromagnet for MP separation
- Deployable NMR for in situ water quality monitoring

Compact-NMR



University of South Carolina



THANKS!



National Science Foundation

Compact-NMR (cNMR)



Our design is open source and available on GitHub!



<https://github.com/ARTS-Laboratory/Compact-NMR>

This work was supported by the National Science Foundation under Grant Nos. 2037771 and 2237696. The authors extend their gratitude to the National Science Foundation for such support. Discussed findings are those of the authors and do not necessarily reflect the views of the National Science Foundation.