

Predicting Sustainable Aviation Fuel Mixtures using Low-Resolution Nuclear Magnetic Resonance

Winford Janvrin¹, Parker Huggins^{1,2}, Austin Downey^{1,3}, Sang Hee Won¹

¹Department of Mechanical Engineering, University of South Carolina

²Department of Electrical Engineering, University of South Carolina

³Department of Civil and Environmental Engineering, University of South Carolina



UNIVERSITY OF
SOUTH CAROLINA

Introduction


Materials &
Methods

Results &
Discussion

Conclusion

Aviation's Environmental Impact


- The Aviation industry accounted for 2.4% of total green house gas emissions in 2018.
- The demand for jet fuel is projected to grow from 71 billion gallons to over 230 billion gallons by 2050, emphasizing the need for sustainable solutions.
- Sustainable Aviation Fuel (SAF) offers a transformative opportunity, reducing greenhouse gas emissions by up to 80% compared to traditional jet fuels.



OurWorldinData.org/co2-and-greenhouse-gas-emissions | CC BY

Sustainable Aviation Fuel (SAF)

- Sustainable aviation fuel (SAF) is a biofuel that can be used in place of conventional jet fuel to power aircraft
- SAF offers a transformative opportunity, reducing greenhouse gas emissions by up to 80% compared to traditional jet fuels.



Svensson, Christian, Amir AM Oliveira, and Tomas Grönstedt. "Hydrogen fuel cell aircraft for the Nordic market." *International Journal of Hydrogen Energy* 61 (2024): 650-663.

Introduction

Materials &
Methods

Results &
Discussion

Conclusion


The Case for SAF Detection

- The United States has launched the Sustainable Aviation Fuel Grand Challenge aiming to accelerate the production of SAF to meet 100% of commercial demand by 2050.
- Accurate detection of SAF ensures:
 - Compliance with stringent regulations.
 - Achievement of intended environmental benefits.
 - Optimization of engine performance.
 - Supporting widespread adoption in the aviation industry.



Rolls-Royce engine compatible with 100% Sustainable Aviation Fuel

Credit: Rolls-Royce




Team Goals (now and future)

1. Rapidly differentiate between mixtures of petroleum-derived and renewable jet fuels.
2. Create a cost-effective on-site fuel monitoring system.
3. Develop sensors for advanced engine control, potentially for dynamic mixtures of SAF and kerosene.


Approach:

- Time domain nuclear magnetic resonance relaxometry (TD-NMR) using a CPMG sequence
- Inexpensive compared to High-field NMR spectroscopy and portable.





Typical Current Engine Control

Fundamentals of Aircraft Turbine Engine Control, Dr. Sanjay Garg Chief, Controls and Dynamics Branch, NASA




ARTS-Lab desktop NMR system

- Control handled by LabVIEW program and NI-PXI chassis
- All electronics (barring one amplifiers) housed on a single PCB
- GUI developed for easy data acquisition and export




Extraction of T2 relaxation curve using a CPMG sequence




Permanent magnet array


- N42 cylindrical dipole magnets enclosed by a steel yolk
- 1018 carbon steel caps affixed to magnet surfaces
- Peak flux density of 0.645 T \rightarrow Larmor frequency of 27.5 MHz
- Temperature shift gradient of -800 ppm/K



Simulation of magnet flux density




Magnet dimensions



RF electronics

- A single 24 V DC power supply required
- Impedance of all cables and PCB traces matched to 50Ω
- Waveform generator \rightarrow sine wave at Larmor frequency
- Pulse generator \rightarrow CPMG pulse train
- Duplexer (crossed diodes) isolates probe and LNA


General flow







Signal generation and control

- NI PXI chassis
 - Arbitrary waveform generator
 - Pulse train generator
 - 16-bit digitizer



- CPMG pulse train
 - 3955 total pulses
 - 90° pulse duration is 7 μs
 - $\tau = 0.625$ ms






Data acquisition

- LabVIEW GUI serves as front end
- Each test comprises 8 scans (average)
- Time for T_2 curve acquisition < 1 min
- Thermocouple used for frequency calibration



current scan outputs



user adjustable parameters

TD-NMR signal and T_2


- T_2 relaxation modeled as $M_{xy}(t) = M_0 \exp(-t/T_2)$
- Relaxation rate is the reciprocal of relaxation time (i.e., $R_2 = 1/T_2$)
- Linear relationship between R_2 and hydrogen content well established



Sample preparation


- Fuel mixture sets of Jet A and HRJ Camelina were created in 10% mass increments, with each mixture totaling 0.3 grams.
 - Ex. 0.27 grams Jet A, 0.03 grams HRJ Camelina
- 11 distinct mixtures including pure samples of Jet A and HRJ Camelina were probed five times, generating a dataset of 55 T_2 curves.






Relaxation rate analysis


- The averaged T_2 data is shown in figure (a) for 6 different concentrations of Jet A and HRJ Camelina.
- It was observed that the relaxation rates decrease with increasing concentrations of Jet A.





SAF Relaxation rate vs concentration

- Relaxation rates decrease with increasing Jet A concentration, reflecting a link between composition and TD-NMR response.
- A strong linear correlation ($R^2 = 0.9845$) between measured and synthetic relaxation rates confirms TD-NMR's reliability for quantifying SAF concentrations.





Predicting mixtures


- Synthetic relaxation curves were generated by superimposing the decay curves of Jet A and HRJ Camelina as:

$$M_{\text{weighted}}(t) = C_{\text{Jet}} \cdot M_{\text{Jet}}(t) + C_{\text{HRJ}} \cdot M_{\text{HRJ}}(t),$$

where $M_{\text{Jet}}(t)$ is the relaxation curve of Jet A, $M_{\text{HRJ}}(t)$ is the relaxation curve of HRJ Camelina, and $0 \leq C_{\text{Jet}}, C_{\text{HRJ}} \leq 1$ are the concentrations of Jet A and HRJ Camelina, respectively.

- The top figure shows each sample's measured R_2 value on the horizontal axis and the corresponding synthesized value on the vertical axis. A linear relationship was found achieving an R^2 of 0.9845.
- Fitting error is shown in the bottom figure and is attributed to human error during preparation and mixing small volumes of fuel.





Conclusion

- Time-domain nuclear magnetic resonance (TD-NMR) is a cost-effective and reliable tool for analyzing Sustainable Aviation Fuel (SAF) mixtures.
- The system developed is compact, open-source, and adaptable, enabling widespread use across industries.
- High accuracy and linear correlation demonstrate the potential for real-time, on-site SAF monitoring.
- Future work includes integrating flow-through systems for continuous monitoring and exploring broader applications like material science and food safety.



Credit: Envato Elements CC

THANKS!



National Science Foundation

Compact-NMR (cNMR)



Our design is open source and available on GitHub!



<https://github.com/ARTS-Laboratory/Compact-NMR>

This material was sponsored by the ARO under Grant No. W911NF21-1-0306 and the National Science Foundation under Grant Nos. 2152896, 2237696, and 2344357. The views and conclusions contained within this document are those of the authors and should not be interpreted as representing the official policies of the ARO, NSF, or U.S. Government.